

Benzoin


Catalog No: #S3939

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Benzoin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	212.24
Formulation	C14H12O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C1=CC=C(C=C1)C(C(=O)C2=CC=CC=C2)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Alcohol: Soluble Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only