

Choline Fenofibrate

Catalog No: #S3941


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Choline Fenofibrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	421.91
Formulation	C22H28ClNO5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: CC(C)(C(=O)[O-])OC1=CC=C(C=C1)C(=O)C2=CC=C(C=C2)Cl.C[N+](C)(C)CC Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 6.2 mg/mL (Need ultrasonic and warming) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Jones PH, et al. Evaluation of a new formulation of fenofibric acid, ABT-335, co-administered with statins : study design and rationale of a phase III clinical programme. Clin Drug Investig. 2008;28(10):625-34.

Note: This product is for in vitro research use only