

LPAAT- δ Polyclonal Antibody

Catalog No: #41110

Package Size: #41110-1 50ul #41110-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LPAAT- δ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	LPAAT- δ Polyclonal Antibody detects endogenous levels of LPAAT- δ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human LPAAT- δ .
Target Name	LPAAT- δ
Other Names	AGPAT4; 1-acyl-sn-glycerol-3-phosphate acyltransferase delta; 1-acylglycerol-3-phosphate O-acyltransferase 4; 1-AGP acyltransferase 4; 1-AGPAT 4; Lysophosphatidic acid acyltransferase delta; LPAAT-delta
Accession No.	Swiss-Prot: Q9NRZ5NCBI Gene ID: 56895
Uniprot	Q9NRZ5
GeneID	56895;
SDS-PAGE MW	44kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/40000.

Not yet tested in other applications.

Images

Western Blot analysis of MCF7 cells using LPAAT- δ Polyclonal Antibody

Note: This product is for in vitro research use only