

Salermide

Catalog No: #S3953


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Salermide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Sirtuin inhibitor
Calculated MW	394.474
Formulation	C ₂₆ H ₂₂ N ₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC(C1=CC=CC=C1)C(=O)NC2=CC=CC(=C2)N/C=C/3C(=O)C=CC4=CC=CC=C4 Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair Receptor: SIRT1; SIRT2 Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Lara E, et al. Salermide, a Sirtuin inhibitor with a strong cancer-specific proapoptotic effect. *Oncogene*. 2009 Feb 12;28(6):781-91.

Note: This product is for in vitro research use only