

## Adelmidrol

Catalog No: #S3954


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Adelmidrol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	274.361
Formulation	C13H26N2O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:C(CCCC(=O)NCCO)CCCC(=O)NCCOPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 34 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Petrosino S,et al. Adelmidrol increases the endogenous concentrations of palmitoylethanolamide in canine keratinocytes and down-regulates an inflammatory reaction in an in vitro model of contact allergic dermatitis. Vet J. 2016 Jan;207:85-91.

Note: This product is for in vitro research use only