


Description

Product Name	WAY 200070
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor inhibitor
Calculated MW	422.573
Formulation	C25H34N4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: cardiovascular

system SMILES: COC1=CC=CC=C1N2CCN(CC2)CCN(C3=CC=CC=N3)C(=O)C4CCCCC4 Pathways: Endocrinology/Hormones Receptor: ERRα;

ERRβ Boiling pt: Melting pt: Solubility: DMSO: 31mg/ML Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Hughes ZA, et al. WAY-200070, a selective agonist of estrogen receptor beta as a potential novel anxiolytic/antidepressant agent. *Neuropharmacology*. 2008 Jun;54(7):1136-42.

Note: This product is for in vitro research use only