


Description

Product Name	TPI-1
Brief Description	Inhibitors
Purification	99.00%
Target Name	Glucocorticoid Receptor inhibitor
Calculated MW	253.078
Formulation	C ₁₂ H ₆ Cl ₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C1=CC(=C(C=C1Cl)C2=CC(=O)C=CC2=O)Cl Pathways: Endocrinology/Hormones Receptor: recombinant SHP-1 Binding pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kundu S, et al. Novel SHP-1 inhibitors tyrosine phosphatase inhibitor-1 and analogs with preclinical anti-tumor activities as tolerated oral agents. *J Immunol.* 2010 Jun 1;184(11):6529-36.

Note: This product is for in vitro research use only