

LPAAT- $\theta$  Polyclonal Antibody

Catalog No: #41111

Package Size: #41111-1 50ul #41111-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LPAAT- $\theta$ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	LPAAT- $\theta$ Polyclonal Antibody detects endogenous levels of LPAAT- $\theta$ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human LPAAT- $\theta$ .
Target Name	LPAAT- $\theta$
Other Names	AGPAT9; GPAT3; MAG1; HMFN0839; Glycerol-3-phosphate acyltransferase 3; GPAT-3; 1-acylglycerol-3-phosphate O-acyltransferase 9; 1-AGP acyltransferase 9; 1-AGPAT 9; Acyl-CoA:glycerol-3-phosphate acyltransferase 3; hGPAT3; Lung cancer metastas
Accession No.	Swiss-Prot: Q53EU6NCBI Gene ID: 84803
Uniprot	Q53EU6
GeneID	84803;
SDS-PAGE MW	48kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Note: This product is for in vitro research use only