

## Casin

Catalog No: #S3971


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Casin
Brief Description	Inhibitors
Purification	98.00%
Target Name	STAT inhibitor
Calculated MW	306.409
Formulation	C <sub>20</sub> H <sub>22</sub> N <sub>2</sub> O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: C<sub>1</sub>CC(C<sub>2</sub>=C(C<sub>1</sub>)C<sub>3</sub>=C(N<sub>2</sub>)C=CC(=C<sub>3</sub>)C<sub>4</sub>=CC=CC=C<sub>4</sub>)NCCO Pathways: JAK/STAT signaling; Stem Cells Receptor: STAT3 Boling pt: Melting pt: Solubility: DMSO: 10 mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Li Y, et al. Suppression of cancer relapse and metastasis by inhibiting cancer stemness. Proc Natl Acad Sci U S A. 2015 Feb 10;112(6):1839-44.

Note: This product is for in vitro research use only