

LRAT Polyclonal Antibody

Catalog No: #41112

Package Size: #41112-1 50ul #41112-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LRAT Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	LRAT Polyclonal Antibody detects endogenous levels of LRAT protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human LRAT.
Target Name	LRAT
Other Names	LRAT; Lecithin retinol acyltransferase; Phosphatidylcholine--retinol O-acyltransferase
Accession No.	Swiss-Prot: O95237NCBI Gene ID: 9227
Uniprot	O95237
GeneID	9227;
SDS-PAGE MW	27kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of COLO205 cells using LRAT Polyclonal Antibody



Note: This product is for in vitro research use only