

Cytosporone B

Catalog No: #S3976


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cytosporone B
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others agonist
Calculated MW	322.401
Formulation	C18H26O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CCCCCCCC(=O)C1=C(C=C(C=C1CC(=O)OCC)O)O Pathways: Others Receptor: Nur77 Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhan, Y. , et al. Cytosporone B is an agonist for nuclear orphan receptor Nur77. Nature Chemical Biology 4, 548-556 (2008).

Note: This product is for in vitro research use only