

Emedastine

Catalog No: #S3979


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Emedastine |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Histamine Receptor antagonist |
| Calculated MW | 302.422 |
| Formulation | C17H26N4O |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Others SMILES: CCOCCN1C2=CC=CC=C2N=C1N3CCCN(CC3)C Pathways: GPCR/G Protein; Immunology/Inflammation; Neuroscience Receptor: H1 receptor Boiling pt: Melting pt: Solubility: Water: 1.44 mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Liu RF., Efficacy of olopatadine hydrochloride 0.1%, emedastine difumarate 0.05%, and loteprednol etabonate 0.5% for Chinese children with seasonal allergic conjunctivitis: a randomized vehicle-controlled study.

Note: This product is for in vitro research use only