

Vulpinic Acid

Catalog No: #S3982


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vulpinic Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	322.316
Formulation	C19H14O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: COC(=O)/C(=C/1C(=C(C(=O)O1)C2=CC=CC=C2)O)/C3=CC=CC=C3 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Koparal AT., et al. Anti-angiogenic and antiproliferative properties of the lichen substances (-)-usnic acid and vulpinic acid. *Z Naturforsch C.* 2015;70(5-6):159-64.

Note: This product is for in vitro research use only