

Verinurad

Catalog No: #S3994


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Verinurad
Brief Description	Inhibitors
Purification	98.00%
Target Name	OAT inhibitor
Calculated MW	348.42
Formulation	C20H16N2O2S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:MetabolismSMILES:CC(C)(C(=O)O)SC1=C(C=NC=C1)C2=CC=C(C3=CC=CC=C3)C#NPathways:Membrane transporter/Ion channelReceptor:URAT1Boiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Diaz-Torn C, Perez-Herrero N, Perez-Ruiz F. New medications in development for the treatment of hyperuricemia of gout. *Curr Opin Rheumatol.* 2015 Mar;27(2):164-9.

Note: This product is for in vitro research use only