

Amcasertib

Catalog No: #S3995

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Amcasertib
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	539.69
Formulation	C31H33N5O2S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research

Area:CancerSMILES:CCN(CC)CCNC(=O)C1=C(NC(=C1C)/C=C\2/C3=C(C=CC(=C3)C4=CSC(=N4)C5=CC=CC=C5)NC2=O)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 30 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only