

Dapiprazole Hydrochloride

Catalog No: #S3996


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dapiprazole Hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor inhibitor
Calculated MW	361.92
Formulation	C19H28ClN5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC1=CC=CC=C1N2CCN(CC2)CCC3=NN=C4N3CCCC4. CIPathways: GPCR/G Protein Receptor: α 1-adrenergic receptor Boiling pt: Melting pt: Solubility: DMSO: 31 mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Hogan TS., et al. Dose-response study of dapiprazole HCl in the reversal of mydriasis induced by 2.5% phenylephrine. J Ocul Pharmacol Ther. 1997 Aug;13(4):297-302

Note: This product is for in vitro research use only