


Description

Product Name	WEHI-345
Brief Description	Inhibitors
Purification	98.00%
Target Name	TNF inhibitor
Calculated MW	401.2
Formulation	C22H23N7O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:O=C(NCC(C)(N1N=C(C2=CC=C(C)C=C2)C3=C(N)N=CN=C31)C)C4=CC=NC=C4Pathways:ApoptosisReceptor:RIPK2Boiling pt:Melting pt:Solubility:DMSO: 32 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Nachbur U, et al. A RIPK2 inhibitor delays NOD signalling events yet prevents inflammatory cytokine production. Nat Commun. 2015 Mar 17; 6:6442.

Note: This product is for in vitro research use only