

()-Equol

Catalog No: #S4000


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	()-Equol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	242.3
Formulation	C15H14O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:OC1=CC=C2CC(C3=CC=C(O)C=C3)COC2=C1Pathways:Endocrinology/HormonesReceptor:ER α ; ER β Boling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Mahalingam S, Gao L, Gonnering M, Helferich W, et al.Toxicol Appl Pharmacol. 2016 Mar 15;295:47-55.

Note: This product is for in vitro research use only