

Pentostatin

Catalog No: #S4006


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pentostatin
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR inhibitor
Calculated MW	268.27
Formulation	C ₁₁ H ₁₆ N ₄ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: C1[C@@H]([C@H](O[C@H]1N2C=NC3=C2NC=NC[C@H]3O)CO)O Pathways: Neuroscience Receptor: Adenosine receptor Boiling pt: Melting pt: Solubility: Water: 100 mM; DMSO: 75 mM Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Nabhan, et Al. Purine nucleoside analogues and combination therapies in B-cell chronic lymphocytic leukemia: dawn of a new era. Leukemia Res. 28 429.

Note: This product is for in vitro research use only