

## LYAR Polyclonal Antibody

Catalog No: #41115

Package Size: #41115-1 50ul #41115-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	LYAR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	LYAR Polyclonal Antibody detects endogenous levels of LYAR protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human LYAR.
Target Name	LYAR
Other Names	LYAR; PNAS-5; Cell growth-regulating nucleolar protein
Accession No.	Swiss-Prot: Q9NX58NCBI Gene ID: 55646
Uniprot	Q9NX58
GeneID	55646;
SDS-PAGE MW	44kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of SH-SY5Y cells using LYAR Polyclonal Antibody



Western Blot analysis of HELA Jurkat cells using LYAR Polyclonal Antibody

Note: This product is for in vitro research use only