

Relebactam

Catalog No: #S4017


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Relebactam
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	348.38
Formulation	C ₁₂ H ₂₀ N ₄ O ₆ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: C1C[C@H](N2C[C@@H]1N(C2=O)OS(=O)(=O)O)C(=O)NC3CCNCC3 Pathways: Others Receptor: diazabicyclooctane Boiling pt: Melting pt: Solubility: 10 mM in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lapuebla A, et al. Activity of Imipenem with Relebactam against Gram-Negative Pathogens from New York City. Antimicrob Agents Chemother. 2015 Aug;59(8):5029-31.

Note: This product is for in vitro research use only