

Cadazolid

Catalog No: #S4020


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cadazolid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	585.19
Formulation	C29H29F2N3O8
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:O=C(C1=CN(C2CC2)C3=C(C=C(F)C(N4CCC(O)(COC5=CC=C(N6C(O[C@@H](CO)C6)=O)C=C5F)CC4)=C3)C1=O)OPathway
s:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C
and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Seiler P., et al. Cadazolid does not promote intestinal colonization of vancomycin-resistant enterococci in mice. Antimicrob Agents Chemother. 2015 Oct 26;60(1):628-31.

Note: This product is for in vitro research use only