

LYPLA1 Polyclonal Antibody

Catalog No: #41116

Package Size: #41116-1 50ul #41116-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LYPLA1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	LYPLA1 Polyclonal Antibody detects endogenous levels of LYPLA1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human LYPLA1.
Target Name	LYPLA1
Other Names	LYPLA1; APT1; LPL1; Acyl-protein thioesterase 1; APT-1; hAPT1; Lysophospholipase 1; Lysophospholipase I; LPL-I; LysoPLA I
Accession No.	Swiss-Prot: O75608NCBI Gene ID: 10434
Uniprot	O75608
GeneID	10434;
SDS-PAGE MW	25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using LYPLA1 Polyclonal Antibody

Note: This product is for in vitro research use only