

L-Argininamide dihydrochloride

Catalog No: #S4023


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	L-Argininamide dihydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	246.14
Formulation	C6H17Cl2N5O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C(C[C@@H](C(=O)N)N)CN=C(N)N.Cl.Cl Pathways: Others Receptor: ligand binding DNA aptamers Boiling pt: Melting pt: Solubility: Water Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Nutiu R., et al. Aptamers with fluorescence-signaling properties. *Methods*. 2005 Sep;37(1):16-25.

Note: This product is for in vitro research use only