

KIN1148

Catalog No: #S4027

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KIN1148
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	361.44
Formulation	C19H11N3OS2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Others SMILES: O=C(NC1=NC2=CC=C3C(N=CS3)=C2S1)C4=CC5=C(C=C4)C=CC=C5 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 100 mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Probst P., et al. A small-molecule IRF3 agonist functions as an influenza vaccine adjuvant by modulating the antiviral immune response.

Note: This product is for in vitro research use only