

## S1p-receptor-agonist-1

Catalog No: #S4031


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	S1p-receptor-agonist-1
Brief Description	Inhibitors
Purification	98.00%
Target Name	S1P Receptor inhibitor
Calculated MW	409.45
Formulation	C23H24FN3O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:CN1CCN(CC1)CC(=O)N(C)C2=CC=C(C=C2)N/C(=C3/C4=C(C=C(C=C4)C(=O)OC)NC3=O)/C5=CC=CC=C5Pathways:GPCR/G  
 ProteinReceptor:S1PBoling pt:Melting pt:Solubility:DMSO: 6.8 mg/mL (Need ultrasonic and warming)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.WO 2015039587 A1

Note: This product is for in vitro research use only