


## Description

Product Name	MSX-127
Brief Description	Inhibitors
Purification	98.00%
Target Name	CXCR inhibitor
Calculated MW	308.38
Formulation	C16H24N2O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:C1COCCN1CC2=CC(=C(C=C2O)CN3CCOCC3)OPathways:Autophagy; GPCR/G Protein;

Immunology/InflammationReceptor:CXCR4Boiling pt:Melting pt:Solubility:DMSO: SolubleAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. CXCR4 antagonists for the treatment of HIV infection

Note: This product is for in vitro research use only