

Talabostat mesylate

Catalog No: #S4042


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Talabostat mesylate
Brief Description	Inhibitors
Purification	98.00%
Target Name	DPP-4 inhibitor
Calculated MW	310.18
Formulation	C10H23BN2O6S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:B([C@@H]1CCCN1C(=O)[C@H](C(C)C)N)(O)O.CS(=O)(=O)OPathways:Proteases/ProteasomeReceptor:DPP-4Boling
 pt:Melting pt:Solubility:DMSO: 40 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Cristillo AD, et al. HIV-1 Env vaccine comprised of electroporated DNA and protein co-administered with Talabostat. Biochem Biophys Res Commun. 2008 May 23;370(1):22-6.

Note: This product is for in vitro research use only