

## Etoxazole

Catalog No: #S4055


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Etoxazole
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	359.41
Formulation	C <sub>21</sub> H <sub>23</sub> F <sub>2</sub> NO <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:CCOC1=C(G=CC(=C1)C(C)(C)C)C2COC(=N2)C3=C(C=CC=C3F)FPathways:OthersReceptor:chitin biosynthesisBoiling pt:Melting pt:Solubility:chloroform; Ethyl acetate; MethanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Uesugi R, et al. Genetic basis of resistances to chlorfenapyr and etoxazole in the two-spotted spider mite (Acari: Tetranychidae). J Econ Entomol. 2002 Dec;95(6):1267-74.

Note: This product is for in vitro research use only