

CYCLIZINE HYDROCHLORIDE

Catalog No: #S4056


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CYCLIZINE HYDROCHLORIDE
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	302.846
Formulation	C18H22N2·HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CN1CCN(CC1)C(C2=CC=CC=C2)C3=CC=CC=C3.Cl Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO (Slightly, Heated), Methanol (Slightly, Heated), Water (Slightly, Sonicated) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kuntzman, R., et al.: J. Pharmacol. Exp. Ther., 158, 332 (1967)

Note: This product is for in vitro research use only