

beta-Cyclodextrin methyl ethers

Catalog No: #S4072


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	beta-Cyclodextrin methyl ethers
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	1310
Formulation	C54H94O35
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:CO[C@H]1[C@@H]([C@H]2[C@H](O[C@@H]1O[C@@H]3[C@H](O[C@@H]([C@H]([C@@H]3OC)OC)O[C@@H]4[C@H](O[C@@H]([C@H]([C@@H]4OC)OC)O[C@@H]5[C@H](O[C@@H]([C@H]([C@@H]5OC)OC)O[C@@H]6[C@H](O[C@@H]([C@H]([C@@H]6OC)OC)O[C@@H]7[C@H](O[C@@H]([C@H]([C@@H]7OC)OC)O[C@@H]8[C@H](O[C@H](O2)[C@H]([C@@H]8O)OC)CO)CO)CO)CO)CO)CO)OPathways:OthersReceptor:OthersBoling pt:Melting pt:Solubility:DMSO: 31 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.NA

Note: This product is for in vitro research use only