

Goserelin acetate

Catalog No: #S4102


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Goserelin acetate
Brief Description	Inhibitors
Purification	98.00%
Target Name	GNRH Receptor inhibitor
Calculated MW	1329.46
Formulation	C??H??N??O??
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:CC(C)C[C@@H](C(=O)N[C@@H](CCCN=C(N)N)C(=O)N1CCC[C@@H]1C(=O)NNC(=O)N)NC(=O)[C@@H](COC(C)(C)C)NC(=O)[C@@H](CC2=CC=C(C=C2)O)NC(=O)[C@@H](CO)NC(=O)[C@@H](CC3=CNC4=CC=CC=C43)NC(=O)[C@@H](CC5=CN=CN5)NC(=O)[C@@H]6CCC(=O)N6.CC(=O)OPathways:GPCR/G ProteinReceptor:GNRHRBoiling pt:Melting pt:Solubility:DMSO: 28 mg/MIAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.1.K. Morton; Judith M. Hall (6 December 2012).?Concise Dictionary of Pharmacological Agents: Properties and Synonyms. Springer Science & Business Media. pp.?136C.?

Note: This product is for in vitro research use only