

Complanatoside

Catalog No: #S4130


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Complanatoside
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	624.55
Formulation	C28H32O16
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:COC1=CC(=C2C(=C1)OC(=C(C2=O)O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O)C4=CC=C(C=C4)O[C@H]5[C@@H]([C@H]([C@@H]([C@H](O5)CO)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 100 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Ng YF, et al. Semen Astragali Complanati: an ethnopharmacological, phytochemical and pharmacological review. J Ethnopharmacol. 2014 Aug 8;155(1):39-53.

Note: This product is for in vitro research use only