

LAPACHOL

Catalog No: #S4136


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	LAPACHOL
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	242.274
Formulation	C15H14O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CC(=CCC1=C(C2=CC=CC=C2C(=O)C1=O)O)CPathways:OthersReceptor:othersBoiling pt:Melting pt:Solubility:ethanol: soluble 10 mg/mL clear, light yellow to yellowAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Costa WF,et al. Lapachol as an epithelial tumor inhibitor agent in Drosophila melanogaster heterozygote for tumor suppressor gene wts.Genet Mol Res. 2011 Dec 22;10(4):3236-45.

Note: This product is for in vitro research use only