

4,7-DIHYDROXYCOUMARIN

Catalog No: #S4137


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	4,7-DIHYDROXYCOUMARIN
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	178.143
Formulation	C9H6O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:C1=CC2=C(C=C1O)OC(=O)C=C2OPathways:OthersReceptor:othersBoiling pt:Melting

pt:Solubility:N/AAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only