

Repertaxin

Catalog No: #S4163


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Repertaxin
Brief Description	Inhibitors
Purification	98.00%
Target Name	CXCR inhibitor
Calculated MW	283.39
Formulation	C14H21NO3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
SMILES: CC1=CC(=C(C(=C1)C)S(=O)(=O)N[C@H](C=O)C(C)C)C
Pathways: Autophagy; GPCR/G Protein;
Immunology/Inflammation Receptor: CXCR1; CXCR2
Boiling pt: Melting pt: Solubility: DMSO: 500 mg/mL
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Moriconi A, et al. Design of noncompetitive interleukin-8 inhibitors acting on CXCR1 and CXCR2. *J Med Chem.* 2007 Aug 23;50(17):3984-4002.

Note: This product is for in vitro research use only