


## Description

Product Name	CY-09
Brief Description	Inhibitors
Purification	98.00%
Target Name	NOD inhibitor
Calculated MW	423.43
Formulation	C??H??F?NO?S?
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:N/APathways:Immunology/Inflammation; NF-κBReceptor:NLRP3Boiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Jiang H, et al. Identification of a selective and direct NLRP3 inhibitor to treat inflammatory disorders. J Exp Med. 2017 Nov 6;214(11):3219-3238.

Note: This product is for in vitro research use only