

## B-Raf inhibitor 1 dihydrochloride

Catalog No: #S4167

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	B-Raf inhibitor 1 dihydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Raf inhibitor
Calculated MW	478.94
Formulation	C <sub>26</sub> H <sub>19</sub> ClN <sub>8</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Product Description

## Research

Area: Others SMILES: CC1=C(C2=C(C=C1)C(=NC=C2)NC3=CC=C(C=C3)Cl)NC4=C(C=CC=N4)C5=C6C(=NC=N5)N=CN6.Cl.Cl Pathways: MAPK Receptor: B-Raf Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Wang X, et al. J Med Chem. 2012 Sep 13;55(17): 7332-41.

Note: This product is for in vitro research use only