

# Probucol disuccinate

Catalog No: #S4179


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Probucol disuccinate
Brief Description	Inhibitors
Purification	98.00%
Target Name	LDL inhibitor
Calculated MW	716.987
Formulation	C39H56O8S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

### Research

Area:OthersSMILES:CC(C)(C)C1=CC(=CC(=C1OC(=O)CCC(=O)O)C(C)(C)C)SC(C)(C)SC2=CC(=C(C(=C2)C(C)(C)C)OC(=O)CCC(=O)O)C(C)(C)C  
 athways:MetabolismReceptor:LDLBoiling pt:Melting pt:Solubility:N/AAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.Yamamoto A,et al.

A unique antilipidemic drug--probucol. J Atheroscler Thromb. 2008 Dec;15(6):304-5. Epub 2008 Dec 11.

Note: This product is for in vitro research use only