

## Adipamidoxime


Catalog No: #S4180

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Adipamidoxime
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	174.204
Formulation	C6H14N4O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:C(C/C(=N/O)/N)CC/C(=N/O)/N

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:N/AAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only