

Cosan 528

Catalog No: #S4181


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cosan 528
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	262.531
Formulation	C9H9BrClNO?
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: O=C(NC1=CC=C(Br)C=C1)CCl Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

N/A

Note: This product is for in vitro research use only