

Betahistine


Catalog No: #S4195

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Betahistine
Brief Description	Inhibitors
Purification	99.87%
Target Name	Histamine Receptor inhibitor
Calculated MW	136.19
Formulation	C8H12N2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtitisSMILES:CNCCC1=NC=CC=C1Pathways:GPCR/G Protein; Immunology/Inflammation; NeuroscienceReceptor:H1 receptorBoiling pt:Melting pt:Solubility:Soluble in DMSO, not in waterAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

- 1.Albu S,et al . Treatment of Meniere's disease with intratympanic dexamethazone plus high dosage of betahistine. Am J Otolaryngol. 2016 May-Jun;37(3):225-30.

Note: This product is for in vitro research use only