

Ledipasvir (acetone)

Catalog No: #S4203


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ledipasvir (acetone)
Brief Description	Inhibitors
Purification	98.00%
Target Name	HCV Protease antagonist
Calculated MW	947.08
Formulation	C??H??F?N?O?
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:inflammationSMILES:CC(C)[C@@H](C(=O)N1CC2(CC2)C[C@H]1C3=NC=C(N3)C4=CC5=C(C=C4)C6=C(C5(F)F)C=C(C=C6)C7=CC8=C(C=C7)N=C(N8)[C@@H]9[C@H]1CC[C@H](C1)N9C(=O)[C@H](C(C)C)NC(=O)OC)NC(=O)OC.CC(=O)CPathways:Microbiology&virology;
Proteases/ProteasomeReceptor:HCVBoiling pt:Melting pt:Solubility:DMSO:100 mg/mL (112.48 mM)Ethanol:100 mg/mL (112.48 mM)Water

References

1. Smith MA, et al. Ann Pharmacother. 2015 Mar;49(3):343-50.

Note: This product is for in vitro research use only