

## Protein degrader 1 hydrochloride

Catalog No: #S4207


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Protein degrader 1 hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	467.025
Formulation	C22H31ClN4O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Others SMILES: CC1=C(SC=N1)C2=CC=C(C=C2)CNC(=O)C3CC(CN3C(=O)C(C(C)C)C)N)O.Cl Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: 10 mM in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Buckley DL, et al. HaloPROTACS: Use of Small Molecule PROTACs to Induce Degradation of HaloTag Fusion Proteins. ACS Chem Biol. 2015 Aug 21;10(8):1831-7.

Note: This product is for in vitro research use only