


Description

Product Name	Targocil
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	475.948
Formulation	C21H22ClN5O4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:CCN(CC)C1=NC2=C(N=NN2C3=CC(=C(C=C3)OC)OC)S(=O)(=O)C4=CC=C(C=C4)ClPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:N/AAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

- 1.Suzuki T, et al. In vitro antimicrobial activity of wall teichoic acid biosynthesis inhibitors against Staphylococcus aureus isolates. Antimicrob Agents Chemother. 2011 Feb;55(2):767-74.

Note: This product is for in vitro research use only