

broxaldine

Catalog No: #S4218


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	broxaldine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antibiotic inhibitor
Calculated MW	421.088
Formulation	C17H11Br2NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: CC1=NC2=C(C=C1)C(=CC(=C2OC(=O)C3=CC=CC=C3)Br)Br Pathways: Microbiology & virology Receptor: antiprotozoal
 Boiling pt: Melting pt: Solubility: DMSO: 42.86 mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only