


## Description

Product Name	Mitapivat
Brief Description	Inhibitors
Purification	98.00%
Target Name	PKM inhibitor
Calculated MW	450.557
Formulation	C24H26N4O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C1CC1CN2CCN(CC2)C(=O)C3=CC=C(C=C3)NS(=O)(=O)C4=CC=CC5=C4N=CC=C5Pathways:MetabolismReceptor:PKM2Boli  
ng pt:Melting pt:Solubility:10 mM in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.C Kung. Activators of pyruvate kinase M2 and methods of treating disease. PCT Int. Appl. WO 2013056153 A1.

Note: This product is for in vitro research use only