

T.cruzi Inhibitor

Catalog No: #S4248

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T.cruzi Inhibitor
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antibiotic inhibitor
Calculated MW	426.91
Formulation	C22H26ClF3N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research

Area: Others SMILES: CN(CC1=CC=CC=C1)C(=O)C2CCCN(C2)CC3=CC=C(C=C3)C(F)(F)F. CIPathways: Microbiology&virology Receptor: Trypanosoma cruzi
Boiling pt: Melting pt: Solubility: DMSO: Soluble Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

N/A

Note: This product is for in vitro research use only