

MEK-1 Polyclonal Antibody

Catalog No: #41133

Package Size: #41133-1 50ul #41133-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MEK-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	MEK-1 Polyclonal Antibody detects endogenous levels of MEK-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human MEK-1 around the non-phosphorylation site of T286.
Target Name	MEK-1
Other Names	MAP2K1; MEK1; PRKMK1; Dual specificity mitogen-activated protein kinase kinase 1; MAP kinase kinase 1; MAPKK 1; MKK1; ERK activator kinase 1; MAPK/ERK kinase 1; MEK 1
Accession No.	Swiss-Prot: Q02750NCBI Gene ID: 5604
Uniprot	Q02750
GeneID	5604;
SDS-PAGE MW	43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of 293T-UV cells using MEK-1 Polyclonal Antibody



Note: This product is for in vitro research use only