

Hispidol

Catalog No: #S4319


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hispidol
Brief Description	Inhibitors
Purification	98.00%
Target Name	TNF inhibitor
Calculated MW	254.24
Formulation	C15H10O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:C1=CC(=CC=C1/C=C\2/C(=O)C3=C(O2)C=C(C=C3)O)OPathways:ApoptosisReceptor:TNF-αBoiling pt:Melting pt:Solubility:DMSO: 100 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Kadayat TM, et al. Discovery and structure-activity relationship studies of 2-benzylidene-2,3-dihydro-1H-inden-1-one and benzofuran-3(2H)-one derivatives as a novel class of potential therapeutics for inflammatory bowel disease. Eur J Med Chem. 2017 Sep 8;137:575-597.

Note: This product is for in vitro research use only